

ALIFAX MICROBIOLOGY LINE



ANTIMICROBIAL SUSCEPTIBILITY TEST


APPLICATION FOR ALIFAX
MICROBIOLOGY INSTRUMENTS

Code	AST EUCAST - CE marked	Gram-negatives			Gram-positives					
		Enterobacteriales	Pseudomonas spp.	Acinetobacter spp.	S. aureus	S. lugdunensis	S. saprophyticus	S. epidermidis	CNS	Enterococcus spp.
1 SI 956-AMK	Amikacin Enterobacteriaceae / EU	•								
2 SI 978-AMK	Amikacin Pseudomonas / EU		•	•						
3 SI 981-AMK	Amikacin Staphylococci / EU				•	•	•	•	•	
4 SI 980-AMC	Amoxicillin-Clavulanic Acid Enterobacteriales / EU (systemic)	•								
5 SI 954-AMP	Ampicillin Enterobacteriaceae / EU	•								
6 SI 955-AMP	Ampicillin Entecococci / EU									•
7 SI 997-AMS	Ampicillin-Sulbactam Enterobacteriaceae / EU	•								
8 SI 957-ATM	Aztreonam Enterobacteriaceae / EU	•								
9 SI 9A1-FEP	Cefepime Enterobacteriales / EU	•								
10 SI 959-CTX	Cefotaxime / EU	•								
11 SI 949-CAZ	Ceftazidime Enterobacteriaceae / EU	•								
12 SI 950-CAZ	Ceftazidime Pseudomonas / EU		•							
13 SI 988-CZA	Ceftazidime-Avibactam Enterobacteriales / EU	•								
14 SI 951-CRO	Ceftriaxone / EU	•								
15 SI 960-FUR	Cefuroxime / EU	•								
16 SI 961-FOX	Cefoxitin Staph. aureus / EU				•	•	•			
17 SI 962-FOX	Cefoxitin CNS / EU							•	•	
18 SI 963-CIP	Ciprofloxacin / EU	•								
19 SI 9A2-CIP	Ciprofloxacin Pseudomonas/ EU		•							
20 SI 964-CLI	Clindamycin Staphylococci / EU				•	•	•	•	•	
21 SI 984-CST	Colistin Enterobacteriaceae / EU	•								
22 SI 983-CST	Colistin Pseudomonas / EU		•	•						
23 SI 965-SXT	Cotrimoxazole Enterobacteriaceae / EU	•								
24 SI 982-SXT	Cotrimoxazole Staphylococci / EU				•	•	•	•	•	
25 SI 986-DAP	Daptomycin Staphylococci / EU				•	•	•	•	•	
26 SI 967-GEN	Gentamicin / EU	•	•	•						
27 SI 968-GEN	Gentamicin Staphylococci / EU				•	•	•	•	•	
28 SI 999-GEN	Gentamicin HLAR / EU									•
29 SI 969-LEV	Levofloxacin / EU	•								
30 SI 9A4-LEV	Levofloxacin Pseudomonas / EU		•							
31 SI 970-LZD	Linezolid / EU				•	•	•	•	•	•
32 SI 971-MEM	Meropenem Enterobacteriaceae / EU	•								
33 SI 979-MEM	Meropenem Pseudomonas / EU		•	•						
34 SI 953-TZP	Piperacillin-Tazobactam Enterobacteriaceae / EU	•								
35 SI 952-TZP	Piperacillin-Tazobactam Pseudomonas / EU		•							
36 SI 996-RIF	Rifampicin / EU				•	•	•	•	•	
37 SI 975-TEI	Teicoplanin S. aureus and Enterococci / EU				•					•
38 SI 976-TEI	Teicoplanin CNS / EU					•	•	•	•	
39 SI 973-VAN	Vancomycin S. aureus / EU				•					
40 SI 974-VAN	Vancomycin CNS and Enterococci / EU					•	•	•	•	•

Not for sale in U.S.

Due to regulatory requirements, the products listed herein are available for sale in any relevant Country territory provided that they have been registered according to the applicable domestic law. The products listed herein may have incurred changes with respect to the last available edition.

www.alifax.com



Code	AST CLSI - CE marked	Gram-negatives				Gram-positives				
		Enterobacteriaceae	P. aeruginosa	Acinetobacter spp.	Stenotrophomonas maltophilia	S. aureus	S. lugdunensis	S. epidermidis	Other Staphylococcus spp.	Enterococcus spp.
1 SI 801-AMK	Amikacin / CL	•	•	•		•	•	•	•	
2 SI 803-AMP	Ampicillin Enterobacteriaceae / CL	•								
3 SI 804-AMP	Ampicillin Entecoccci / CL									•
4 SI 805-ATM	Aztreonam Enterobacteriaceae / CL	•								
5 SI 806-ATM	Aztreonam Pseudomonas / CL		•							
6 SI 807-CTX	Cefotaxime Enterobacteriaceae / CL	•								
7 SI 811-CAZ	Ceftazidime Enterobacteriaceae / CL	•								
8 SI 812-CAZ	Ceftazidime Pseudomonas / CL		•	•	•					
9 SI 813-CRO	Ceftriaxone Enterobacteriaceae / CL	•								
10 SI 815-FUR	Cefuroxime Enterobacteriaceae / CL	•								
11 SI 809-FOX	Cefoxitin Staph. aureus / CL					•	•			
12 SI 810-FOX	Cefoxitin CNS / EU							•	•	
13 SI 817-CIP	Ciprofloxacin / CL	•								
14 SI 818-CLI	Clindamycin Staphylococci / CL					•	•	•	•	
15 SI 848-CST	Colistin Enterobacteriaceae / CL	•								
16 SI 819-CST	Colistin Pseudomonas / CL		•	•						
17 SI 820-SXT	Cotrimoxazole / CL	•				•	•	•	•	
18 SI 837-DAP	Daptomycin Staphylococci / CL					•	•	•	•	
19 SI 821-FOS	Fosfomycin / CL	•								•
20 SI 822-GEN	Gentamicin / CL	•	•	•		•	•	•	•	
21 SI 823-LEV	Levofloxacin / CL	•	•	•	•	•	•	•	•	•
22 SI 824-LZD	Linezolid / CL					•	•	•	•	•
23 SI 825-MEM	Meropenem Enterobacteriaceae / CL	•								
24 SI 826-MEM	Meropenem Pseudomonas / CL		•	•						
25 SI 829-OXA	Oxacillin Staph. aureus / CL					•	•			
26 SI 830-TZP	Piperacillin-Tazobactam Enterobacteriaceae / CL	•								
27 SI 831-TZP	Piperacillin-Tazobactam Pseudomonas / CL		•	•						
28 SI 832-TEI	Teicoplanin Staphylococci / CL					•	•	•	•	
29 SI 833-TEI	Teicoplanin Enterococci / CL									•
30 SI 834-VAN	Vancomycin Staph. aureus / CL					•				
31 SI 835-VAN	Vancomycin CNS and Enterococci / CL						•	•	•	•

30 TESTS PER KIT

Each Alifax AST antibiotic kit contains:

- 5 vials of lyophilized antibiotic
- 2 vials of regenerating solution
- 1 Mic CARD



SHELF LIFE

From production: 8-12 months
From reconstitution: 3-7 days

STORAGE CONDITIONS

Refrigerated temperature (+4° / +8° C)

ALIFAX S.r.l.

Via Petrarca, 2/1 - Polverara (PD) - Italy

(+39) 049 0992000 | info@alifax.com | www.alifax.com

Entered in Padova companies register at the n. 04337640280. Company stock € 10.000.000 entirely deposited.

Company with single partner and subjected to direction and coordination of Alifax Holding.